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Referencer	<p>1. Henriksen LO, Faber NR, Moller MF, Nexø E, Hansen AB. Stability of 35 biochemical and immunological routine tests after 10 hours storage and transport of human whole blood at 21°C. Scandinavian Journal of Clinical and Laboratory Investigation. 2014;74(7):603-10.</p> <p>2. Kit insert; HAPT2, 03005593322 v14.0 (c501)</p> <p>3. Kit insert; HAPT2, 05992001190 v8.0 (c702)</p> <p>4. Guder WG, Fonseca-Wollheim Fd, Heil W, Schmitt Y. Quality of Diagnostic Samples, Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine. German Society for Clinical Chemistry and Laboratory Medicine 3rd completely revised edition 2010.</p> <p>5. Hornung N, Kristensen JK. Haptoglobin [updated 21.07.2023; cited 29.08.2025. Available from:  <a href="https://www.sundhed.dk/sundhedsfaglig/laegehaandbogen/undersogelser-og-proever/klinisk-biokemi/blodproever/haptoglobin/">https://www.sundhed.dk/sundhedsfaglig/laegehaandbogen/undersogelser-og-proever/klinisk-biokemi/blodproever/haptoglobin/</a>.</p> <p>6. Barcellini W. Diagnosis of hemolytic anemia in adults [updated 01.08.2025; cited 29.08.2025. Available from:  <a href="https://www.uptodate.com/Haptoglobin_adults">https://www.uptodate.com/Haptoglobin_adults</a>.</p> <p>7. Debraun MR. Overview of hemolytic anemias in children [updated 07.08.2024. Available from:  <a href="https://www.uptodate.com/Haptoglobin_children">https://www.uptodate.com/Haptoglobin_children</a>.</p> <p>8. EFLM Biological Variation Database: Haptoglobin [updated 13.08.2025; cited 29.08.2025. Available from:  <a href="https://biologicalvariation.eu/api/search/by_analyte?format=html&amp;query=Haptoglobin">https://biologicalvariation.eu/api/search/by_analyte?format=html&amp;query=Haptoglobin</a>.</p>
------------	---

 Version	Godkendt	Revisionsinformation
1	22.10.2025	Nyoprettet